



LEIMKUEHLER
ORTHOTIC & PROSTHETIC
CENTER, INC.

Specialists Since 1948

Amherst Office

205 N Leavitt Road
Amherst, OH 44001
P: 440-988-5770
F: 440-988-5776

Sandusky Office

1807 W. Perkins Ave.
Sandusky, OH 44870
P: 419-625-5373
F: 419-625-5377

N. Olmsted Office

28350 Lorain Road
North Olmsted, OH 44070
P: 440-979-9226
F: 440-979-9326

www.leimop.com

Gel Liners

- Gel liners should be changed daily. Alternate between the two you received.
- Insert liner with gel side facing out and wash in the sink with your normal soap.
- After washing, either spray or rub with alcohol to disinfect.
- Then turn the liner with fabric side out to dry
- If you get an allergic reaction to a liner, stop wearing it and call the office.
- The liner can cause a rash and odor if they are not kept clean enough.
- You may experience perspiration at first but, the body will slowly adapt and it will become less over time. If this continues however antiperspirant can be used.
- Do not use excessive skin lotion. It can contaminate the liners, use lotion at night (if desired) after you remove the liner for the day.